## MONTANA Arnica cream

Karel Hadek

Code: 2001

The MONTANA cream, as the name itself suggests already, makes use of the active ingredients present in arnica. Arnica montana is familiar for its ability to stimulate blood circulation in the skin, which may in turn improve lymphatic circulation in the place of application. This way it helps skin detoxification and cell regeneration.

Thank to arnica, this preparation may promote minor injury healing and help in reducing swelling. It soothes the skin, contributes to improving skin regeneration and restores its natural functions. The skin becomes soft and supple after application. Also suitable for treating minor surface injuries. It can be described as a special preparation for combination skin.

The preparation has also proved very effective when treating dry elbows. But it can also be used with success in foot care.

Main ingredients: Vegetable oils: sesame, walnut, wheat germ, etc. Arnica macerate. Lecithin. Beeswax. Emulsifiers. Essential oils: juniper berry, rosemary. Vitamins A, E, F. Beta-carotene, panthenol. Allantoin. Glycerin. Distilled water.

**Application:** Wash your face and apply the required amount of the regenerative cream on the skin, working in gently. The skin is not shiny after application.

INCI: Aqua, Glycine soja Oil, Cetearyl alcohol, Juglans regia Seed Oil, Sesamum indicum Seed Oil, Triticum vulgare Germ Oil, Panthenol, Arnica Montana Flower Extract, Butyrospermum parkii Butter, Glycerin, Sodium cetearyl sulfate, Cera flava, Urea, Tocopheryl acetate, Lactic acid, Lecithin, Retinyl palmitate, Ascorbyl palmitate, Melaleuca alternifolia Leaf Oil, Rosmarinus officinalis Leaf Oil, Juniperus communis Fruit Oil, Allantoin, Thymus vulgaris Oil, Thymus serpyllum Oil, Salvia officinalis Oil, Lavandula hybrida Oil, Eugenia caryophyllus Flower Oil, Citrus limon Peel Oil, Thymol, Beta-carotene

Storage: The stated expiry date can be substantially extended by keeping the product in a dark and cold place (ideally in a fridge at 6–10 °C).

Best before the date indicated on the packaging. Batch number and volume are indicated on the packaging.